

PANTHERA BAMBINO HYBRID



THE PANTHERA BAMBINO HYBRID is a light everyday wheelchair for young children. The Bambino hybrid's lightweight frame, and cambered rear wheels make it an easy to drive active chair for children.

Ideal for children ranging in ages 4 – 8 years old who need a light weight easy to handle chair.

Now also available with abduction front (available only in red frame). Standard chassis, available in 3 widths and two colours (red and black) with custom colour options.

STANDARD FEATURES

- Footrest adjustable in height, angle & depth
- 3" front caster
- 18" rear wheels
- Brakes
- Balancing adjustable axle
- Anti-tip device
- Light weight frame: starting at 10.758 lb (4.89 kg)
- Cambered rear axle

ACCESSORIES (Sold Separately)



Adjustable Push Handle



Side Guards



Push Bar



Cushion

PANTHERA BAMBINO HYBRID SPECIFICATIONS

Seat width	9.5" (24 cm)	10.63" (27 cm)	11.8" (30 cm)
Backrest height	6.7 - 9.8" (17 - 25 cm) , 10.6 - 13.8" (27 - 35 cm)		
Seat height (rear)	13.9" (35 cm)		
Seat height (front)	15" (39 cm)		
Seat angle	9-11.8" (23 -30 cm) adjustable		
Seat depth	4"		
Total length	24.75" (63 cm)		
Total height	22.4 - 25.6" (57 - 65 cm)		
Total width	17.3" (44 cm)	18.5" (47 cm)	19.7" (50 cm)
Max User weight	132 lb (60 kg)		
Total weight	15.27 lb (6.94 kg)	15.77 lb (7.16 kg)	16.28 lb (7.4 kg)
Transfer weight	10.77 lb (4.89 kg)	11.31lb (5.14 kg)	11.81 lb (5.37 kg)



Distributed exclusively in North America by



Triumph Mobility Inc.

www.triumphmobility.com | info@triumphmobility.com | 1.855.546.0711
R_1016_BambinoHybrid-SalesSheet_EN | Printed in Canada