









Adjusted by:









Tel:

E-mail:

Adjusted for:

Date:

<p><b>Colour:</b></p>	<input type="checkbox"/> <b>Red</b> <input type="checkbox"/> <b>Black</b>		<p><b>NOTE!</b> Abduction front and seat width 24 cm only in red!</p>
<p><b>Seat width (BB):</b></p> 	<input type="checkbox"/> <b>24 cm *</b> <input type="checkbox"/> <b>27 cm</b> <input type="checkbox"/> <b>30 cm</b> <input type="checkbox"/> <b>33 cm</b> 356 24 00G      356 27 00(BL)G      356 30 00(BL)G      356 33 00(BL)G <input type="checkbox"/> <b>24 cm Abd.*</b> <input type="checkbox"/> <b>27 cm Abd.*</b> <input type="checkbox"/> <b>30 cm Abd.*</b> <input type="checkbox"/> <b>33 cm Abd.*</b> 357 24 00G      357 27 00G      357 30 00G      357 33 00G		<p><b>NOTE!</b> Use "BLG" after art.no. when choosing black colour</p> <p><b>* Only available in red</b></p> <p>Width 24 cm has the back moved 5 cm forward      Width 24-33 cm <b>max. user weight 60 kg</b>                  art.no. includes: seat, backrestlock, footrest, calf strap, brake and anti tip device.</p>
<p><b>Rear axle:</b></p> 	<input type="checkbox"/> <b>Camber 6°</b> 336 BB 16      BB is replaced by seat width, see "Seat width, above" <input type="checkbox"/> <b>other rear axle</b> <a href="#">See list of options.</a> <input type="text"/>		<p><b>Rear axle adjustment</b> balance of the chair</p> <input type="checkbox"/> <b>standard position, 25 mm</b> <input type="checkbox"/> <b>other position</b> Note measure (A) for desired position. A= _____ mm 
<p><b>Castors:</b></p> 	<input type="checkbox"/> <b>5" for 20" rear wheels</b> <input type="checkbox"/> <b>5" thick rubber for 22-24" rear wheels</b> 230 01 92      230 00 96 <input type="checkbox"/> <b>5" for 22-24" rear wheels</b> 230 00 02		<p><b>*NOTE!</b> You need two castors for each chair, the art.nr. is for one castor!</p>
<p><b>Backrest height:</b></p> 	<input type="checkbox"/> <b>Adjustable height 17-25 cm</b> <input type="checkbox"/> <b>Adjustable height 27-35 cm</b> 323 BB 17      323 BB 18		
<p><b>Footrest:</b></p> 	<input type="checkbox"/> <b>Footrest, Bambino</b> <input type="checkbox"/> <b>Footrest, Bambino Swing up</b> (400 21 BB) included in art.nr. 35X BB 00G, see "Seat width"      400 21 10 (assembly-kit for swing-up footrest)		
<p><b>Rear wheels:</b></p> 	<input type="checkbox"/> <b>20", Titanium pushrim</b> <input type="checkbox"/> <b>20", Maxgrip pushrim, thin*</b> 132 00 20      132 00 30 <input type="checkbox"/> <b>22", Titanium pushrim</b> <input type="checkbox"/> <b>22", Maxgrepp, pushrim, thin*</b> 131 00 20      131 00 30 <input type="checkbox"/> <b>24", Titanium pushrim</b> <input type="checkbox"/> <b>24", Spider wheel, friction pushrim</b> 130 00 20      110 00 30 <input type="checkbox"/> <b>24", Tetra grip pushrim</b> <input type="checkbox"/> <b>24", Para grip pushrim</b> 130 00 40      130 00 50 <input type="checkbox"/> <b>24", Spider wheel, titanium pushrim</b> <input type="checkbox"/> <b>24", Maxgrip pushrim, thin*</b> 110 00 20      130 00 30 <input type="checkbox"/> <b>Other rear wheels</b> <a href="#">See list of options.</a> <input type="text"/>		<p><b>*NOTE!</b> Contains latex!</p> <p><b>NOTE!</b> You need two rear wheels for each chair, the art.no. is for one rear wheel!</p>
<p><b>Spoke protectors:</b></p> 	<input type="checkbox"/> <b>24" Spoke protectors</b> <input type="checkbox"/> <b>22" Spoke protectors</b> <input type="checkbox"/> <b>20" Spoke protectors</b> 450 00 00 with logo      450 01 00 with logo      450 00 21 with logo <input type="checkbox"/> <b>24" Spoke protectors</b> <input type="checkbox"/> <b>22" Spoke protectors</b> 450 00 11 without logo      450 01 11 without logo		

<b>Brakes:</b> 	<input type="checkbox"/> <b>High brake, moved ahead</b> (460 00 50) included in art.nr.356 BB 00G, see "Seat width"	sid 2/2
<b>Push handles:</b> 	<input type="checkbox"/> <b>Pushbar</b> max height 100 cm with 22" wheels 422 00 BB max. user weight 60 kg <input type="checkbox"/> <b>Height adjustable +10 cm</b> 425 10 00 max. user weigh 60 kg <input type="checkbox"/> <b>Height adjustable</b> 425 00 00 max. user weigh 60 kg <input type="checkbox"/> <b>Height adjustable +20 cm</b> 425 20 00 max. user weigh 60 kg	
<b>Anti tip device:</b> 	<input type="checkbox"/> <b>Anti tip device, Bambino</b> (435 10 01) included in art.no. 35X BB 00G, see "Seat width"	
<b>Sideguards:</b> 	<input type="checkbox"/> <b>Sideguards Bambino 20"</b> 441 00 15 <input type="checkbox"/> <b>Sideguards Bambino 22"</b> 441 00 10 <input type="checkbox"/> <b>Sideguards Bambino 24"</b> 441 04 10 <input type="checkbox"/> <b>Sideguards Bambino winter 20"</b> 441 01 15 <input type="checkbox"/> <b>Sideguards Bambino winter 22"</b> 441 01 10 <input type="checkbox"/> <b>Sideguards Bambino winter 24"</b> 441 03 10	
<b>Cushions:</b> 	<input type="checkbox"/> <b>Cushion Micro long 2,5 cm for seat width 24 cm</b> 472 41 27 <input type="checkbox"/> <b>Cushion Bambino 2,5 cm</b> 472 30 BB adjustable seatdepth 23-30 cm <input type="checkbox"/> <b>Cushion Micro long 5 cm for seat width 24 cm</b> 475 41 27 <input type="checkbox"/> <b>Cushion Bambino 5 cm</b> 475 30 BB adjustable seatdepth 23-30 cm	
<b>Accessories:</b>	<input type="checkbox"/> <b>Torso strap Bambino</b> 438 02 38 <input type="checkbox"/> <b>Tube protection</b> _____ pcs 990 02 51 protects tube and paint <input type="checkbox"/> <b>Head cushion</b> 471 00 BB	
<b>Back rest angle:</b> adjustment	<b>Back rest angle adjusted with screw</b> Note measure (A) for desired angle      A= _____ mm 	
<b>Footrest:</b> adjustment	<b>Footrest adjustments</b> Specify measures (A and B) for the footrest      A= _____      B= _____ mm 	
<b>Marking:</b>	<input type="checkbox"/> <b>YES</b> mark the wheelchair for transportation in disabled service vehicles. <input type="checkbox"/> <b>NO</b> do <b>NOT</b> mark the wheelchair for transportation in disabled service vehicles <div style="border: 1px solid black; padding: 5px; width: fit-content; margin: 10px auto;"> <b>NB! See instruction for use.</b> </div> 	
<b>Other:</b>		