









Adjusted by:










Tel:

E-mail:

Adjusted for:

Date:

<p>Colour:</p>	<input type="checkbox"/> Black					
<p>Seat width (BB):</p> 	<input type="checkbox"/> 36 cm 347 36 00G	<input type="checkbox"/> 39 cm 347 39 00G	<input type="checkbox"/> 42 cm 347 42 00G	<input type="checkbox"/> 45 cm 347 45 00G <div style="border: 1px solid black; padding: 2px; width: fit-content;"> BB= Seat width, HH= Back rest height </div> <div style="border: 1px solid black; padding: 2px; width: fit-content; margin-top: 5px;"> Art.no includes: seat upholstery, backrestlock </div> <div style="border: 1px solid black; padding: 2px; width: fit-content; margin-top: 5px; float: right;"> Width 36-42 cm max. user weight 100 kg Width 45-50 cm max. user weight 150 kg </div>		
<p>Rear axle:</p> 	<input type="checkbox"/> Carbon, +0 cm 330 23 BB <input type="checkbox"/> Carbon, +1,5 cm 330 23 BC* <input type="checkbox"/> Other rear axle <input type="text"/>	<div style="border: 1px solid black; padding: 2px; width: fit-content;"> The wheelchair as narrow as possible </div> <div style="border: 1px solid black; padding: 2px; width: fit-content; margin-top: 5px;"> When choosing option armrest and / or Spinerig-wheels </div> <div style="border: 1px solid black; padding: 2px; width: fit-content; margin-top: 5px;"> See list of options. </div> <div style="border: 1px solid black; padding: 2px; width: fit-content; margin-top: 5px; float: right;"> * sittbredd BB +1 </div>	<div style="border: 1px solid black; padding: 2px; width: fit-content;"> Rear axle adjustment balance of the chair <input type="checkbox"/> standard position, 25 mm <input type="checkbox"/> Other position A= _____ mm </div> <div style="border: 1px solid black; padding: 2px; width: fit-content; margin-top: 5px; float: right;"> Ange måttet (A) för önskad balansering. </div> 			
<p>Castors:</p> 	<input type="checkbox"/> 5" thick rubber 230 00 43	<input type="checkbox"/> 6" 200 45 15 <div style="border: 1px solid black; padding: 2px; width: fit-content; margin-left: auto;"> NOTE! You need two castors for each chair, the art.no. is for one castor! </div>				
<p>Back rest height:</p> 	<input type="checkbox"/> 20 cm 323 BB 20	<input type="checkbox"/> 25 cm 323 BB 25	<input type="checkbox"/> 30 cm 323 BB 30	<input type="checkbox"/> 35 cm 323 BB 35	<input type="checkbox"/> 40 cm 323 BB 40	<input type="checkbox"/> 45 cm 323 BB 45
<p>Special back rests:</p> 	<input type="checkbox"/> Tapered top available in height 30-40 cm 325 BB HH <input type="checkbox"/> Other special back rests <input type="text"/>	<input type="checkbox"/> Curved back available in height 20-40 cm 324 BB HH <div style="border: 1px solid black; padding: 2px; width: fit-content; margin-left: auto;"> See list of options. </div>				
<p>Footrest:</p> 	<input type="checkbox"/> Footrest foldable 33/36 660 36 00 heel strap not included <input type="checkbox"/> Heel strap 33/36 370 03 35	<input type="checkbox"/> Footrest foldable 39/42 660 39 00 heel strap not included <input type="checkbox"/> Heel strap 39/42 370 03 41	<input type="checkbox"/> Footrest foldable 45 660 45 00 heel strap not included <input type="checkbox"/> Heel strap 45 370 03 47			
<p>Rear wheels:</p> 	<input type="checkbox"/> 24", with Titanium-pushrim 130 00 20 <input type="checkbox"/> 24", with Tetragrip-pushrim 130 00 40 <input type="checkbox"/> Other rear wheels <input type="text"/>	<input type="checkbox"/> 24", with Aluminum-pushrim 130 00 10 <input type="checkbox"/> 24", with Paragrip-pushrim 130 00 50 <input type="checkbox"/> 24", with Maxgrip-pushrim, thin** 110 00 01 <div style="border: 1px solid black; padding: 2px; width: fit-content; margin-left: auto;"> **NOTE! Contains latex! </div> <div style="border: 1px solid black; padding: 2px; width: fit-content; margin-top: 5px;"> NOTE! You need two rear wheels for each chair, the art.no. is for one rear wheel! </div>				
<p>Spoke protectors:</p> 	<input type="checkbox"/> 24" Spoke protectors with logo 450 00 00	<input type="checkbox"/> 24" Spoke protectors without logo 450 00 11				

Brake: 	<input type="checkbox"/> High brake 461 00 00
Push handles: 	<input type="checkbox"/> Fixed 421 00 00 <input type="checkbox"/> Foldable 423 00 00 <input type="checkbox"/> Height adjustable 425 00 00
Anti tip device: 	<input type="checkbox"/> Anti tip device 435 10 00
Sideguards: 	<input type="checkbox"/> Sideguards with cover 442 00 10 <input type="checkbox"/> Sideguards with cover, 26" 442 00 40
Armrest: 	<input type="checkbox"/> Short cushion 23 cm 443 01 00 <input type="checkbox"/> Long cushion 30 cm 443 01 30 <input type="checkbox"/> Extended 35 mm, Short cushion 23 cm 443 01 20 <input type="checkbox"/> Extended 35 mm, Long cushion 30 443 01 40 <div style="border: 1px solid red; padding: 2px; width: fit-content; margin-top: 5px;"> If you choose armrest or / and Spinerger rear wheels you need to choose rear axle Carbon +1,5 cm to increase the space between the wheel and the chassis </div>
Cushion: 	<input type="checkbox"/> Cushion 2,5 cm seatdepth 40cm 472 00 BB <input type="checkbox"/> Cushion 5 cm seatdepth 40cm 475 00 BB <input type="checkbox"/> Cushion 7 cm seatdepth 40cm 477 00 BB
Accessories:	<input type="checkbox"/> Seat bag storage of wallet, keys etc.. 990 02 00 is mounted under the seat <input type="checkbox"/> Crutch holder _____ pcs 990 03 00 NOTE! 2pcs / crutch <input type="checkbox"/> End plug, steel _____ pcs 903 00 11 NOTE! 2pcs / chassis
Back rest angle: adjustment	Back rest angle adjusted with screw Note measure (A) for desired angle. A= _____ mm 
Footrest: adjustment	Footrest adjustments Specify measures (A and B) A= _____ mm B= _____ mm 
Marking:	<input type="checkbox"/> YES mark the wheelchair for transportation in disabled service vehicles <input type="checkbox"/> NO do NOT mark the wheelchair for transportation in disabled service vehicles. <div style="border: 1px solid red; padding: 2px; width: fit-content; margin: 5px auto;"> NB! See Instruction for use. </div> 
Other:	<div style="border: 1px solid black; height: 100px; width: 100%;"></div>