

PANTHERA X[®]



The lightest and most high-tech wheelchair.



THE PANTHERA X is by far the lightest and most high-tech wheelchair ever made. This chair is built from the most sophisticated materials using the most advanced technology available today.

The result is a chair with extremely low weight, but still without compromising stability and comfort. As users ourselves, we can definitely say that **PANTHERA XTM** is a dream come true for active users.

A rigid frame is the key to easy handling. All flex requires energy that comes from the user. By using stiff materials and a box frame design, **PANTHERA XTM** reduces energy losses to nearly zero. **PANTHERA XTM** handles and rolls so effortlessly, you just have to drive it to believe it.

On **PANTHERA XTM**, you can adjust the backrest angle as well as fine-tune the backrest upholstery for optimal comfort and proper back support.

STANDARD FEATURES

- Carbon fibre chassis
- Carbon fibre rear axle
- 24" X rear wheel
- 4 position balancing holes
- Light weight frame: Starting at 4.62 lb (2.1 kg)

PANTHERA X

ACCESSORIES

PANTHERA X™ offers a great variety of high-quality accessories that help enhance the wheelchair comfort and performance. These can be purchased separately,



Carbon Fibre Side Guards

Help protect the user's thighs from coming into contact with the wheels. They can be used to sit on during transfers or as the bump up to reach higher places such as a cupboard. Available in 3 sizes.

Cushion X

Cushion X. Available in width 33 - 45 cm and thickness 5 cm.



Transport Protection

Protects the wheelchair chassis from scrapes during transportation.

Arm Hook

Helps with added stability during extraneous movement.

Tube Protection

Attached around the tubes on the chassis with velcro to prevent damages like scratches and minor dents.

Freewheel

Helps with added stability during navigation on uneven surfaces such as lawns, gravel, and snow.

PANTHERA X SPECIFICATIONS

Seat width	13" (33 cm)	14" (36 cm)	15" (39 cm)	16.5" (42 cm)	17.7" (45 cm)
Backrest height	8.6 - 13.8" (22 - 35 cm)				
Seat height (rear)	17" (43 cm)				
Seat height (front)	18.5" (47 cm)				
Seat angle	7°				
Seat depth	13.8 - 18.1" (35 - 46 cm)				
Total length	32.5" (82.5 cm)				
Total height	25" - 33" (64 - 84 cm)				
Total width	21.25" (54 cm)	22.4" (57 cm)	23.6" (60 cm)	24.8" (63 cm)	26" (66 cm)
Max User weight	220 lb (100 kg)	220 lb (100 kg)	220 lb (100 kg)	220 lb (100 kg)	220 lb (100 kg)
Total weight	9.68 lb (4.4 kg)	9.79 lb (4.45 kg)	9.9 lb (4.5 kg)	10.01 lb (4.55 kg)	10.12 lb (4.6 kg)
Transfer weight	4.62 lb (2.1 kg)	4.73 lb (2.15 kg)	4.84 lb (2.2 kg)	4.95 lb (2.25 kg)	5.06 lb (2.3 kg)



Distributed exclusively
in North America by



Triumph Mobility Inc.

www.triumphmobility.com | info@triumphmobility.com | 1.855.546.0711
R_1116_X-SalesSheet_EN | Printed in Canada