

PANTHERA S3 LARGE



THE PANTHERA S3 LARGE is an active wheelchair designed to let you live life as actively as possible. It is characterised by the attention that we have put into all of its details.

Its low weight combined with the stable fixed frame and casters with a unique profile on the tread surface make this wheelchair very easy to maneuver.

The chassis' seating angle, along with adjustable backrest and seat upholstery and the backrest tilt function combine to provide extremely good seating comfort.

The fabric used in the backrest upholstery has also been carefully selected for ergonomic and comfortable sitting.


The chair's balance point can be adjusted to enable you to find a setting that offers a feeling of security. Accessories such as armrests and side guards have a unique design and you can adjust the height and fore-aft position of the armrests yourself without using tools.

STANDARD FEATURES

- Chrome molybdenum chassis
- Titanium footbar
- 24" Standard rear wheels
- Balancing, continuously adjustable
- Carbon fibre rear axle
- Unique profiled front caster for improved handling.
- Light weight frame: starting at 11.42 lb (5.19 kg)

COLOURS

PANTHERA S3 SHORT Comes standard in a beautiful black frame. It is also available in a series of Premium Colours which can be ordered at an extra cost.



Versatile White



Blue



Mean Green



Panorange



Turtle Turquoise



Brave Bronze

PANTHERA S3 LARGE

ACCESSORIES


PANTHERA S3 LARGE offers a great variety of high-quality accessories that help enhance the wheelchair comfort and performance. These can be purchased separately,



Side Guards with a Soft Upper Section which can be used to sit on during transfers or as the bump up to reach higher places such as a cupboard.



Height-Adjustable Armrest
Adjusts to 4 heights and 4 horizontal positions.



Anti-tip Device for added safety and stability. It can be folded in and out while seated.



Folding Footrest



Height-Adjustable Push Handles

PANTHERA S3 LARGE SPECIFICATIONS

Seat width	15" (39 cm)	16.5" (42 cm)	18" (45 cm)	20" (50 cm)
Backrest height	8" - 18" (20 - 45 cm)			
Seat height (rear)	18" (45.5 cm)			
Seat height (front)	19.5" (49.5 cm)			
Seat angle	7°			
Seat depth	17.7" (45 cm)			
Total length	30.7" -35.4" (78 - 90 cm)			
Total height	26" - 34" (66.5-86.5 cm)			
Total width	23.6" (60 cm)	24.8" (63 cm)	26" (66 cm)	28" (71 cm)
Max User weight	330 lb (150 kg)			
Total weight	20.17 lb (9.17 kg)	20.30 lb (9.23 kg)	20.83 lb (9.46 kg)	20.96 lb (9.53 kg)
Transfer weight	12.35 lb (5.61 kg)	12.48 lb (5.67 kg)	13.00 lb (5.91 kg)	13.13 lb (5.97 kg)



Distributed exclusively
in North America by



Triumph Mobility Inc.

www.triumphmobility.com | info@triumphmobility.com | 1.855.546.0711
R_0617_S3Large-SalesSheet_EN | Printed in Canada