

**ESCAPE Mini Serial Numbers : 2018**

20203- ...

20231- ...

#	Item#	Description	UoM	#	Item#	Description	UoM
1	500-1405	Brake Handle - Complete	EA	34	500-1071	Seat Mechanism - Complete (Assembled)	EA
2	500-1003	Backstrap - Complete	EA	35	500-1090	Soft Seat - LEFT Side	EA
3	500-1004	Backstrap screw & Nut - Complete - 2pc/SET	SET	36	500-1092	Soft Seat - CENTRE	EA
4	500-1002	Backrest connector - Left	EA	37	500-1091	Soft Seat - RIGHT Side	EA
4	500-1005	Backrest Connector - Right	EA	38	500-1095	Seat Handle Strap	EA
6	500-1007	Handle tubing	EA	39	500-1085	Seat Lock Assembly	EA
7	500-1008	Tapping screw (M4*14mm) - 2pc/SET	SET	40	500-1100	Basket hook - RIGHT	EA
8	500-1009	Anatomical handles (1 PR)	PR	41	500-1101	Basket hook - LEFT	EA
9	500-1407	Brake handle lever	EA	42	500-1104	Red clips for shopping bag - LEFT & RIGHT	SET
11	500-1012	Housing / handle bar screw & nut	EA	43	500-1103	Seat to frame screw	EA
13	500-1409	Actuator	EA	44	500-1102	Basket hook screw	EA
14	500-1414	Brake Cable & cover	EA	45	500-1105	Upper folding bar screw	EA
15	500-1415	Handle bar tube Numbered (6 slots) - RIGHT	EA	46	500-1110	Cross bar complete - Escape Super Low	EA
15	500-1416	Handle bar tube Numbered (6 slots) - LEFT	EA	46	500-1115	Cross bar complete - Escape Low	EA
15a	500-1417	Label for handle bar tube - 2pc/SET	SET	46	500-1120	Cross bar complete - Escape Standard	EA
16	500-1425	Adj Collar complete - RIGHT	EA	47	500-1125	Lower cross bar bracket - Super Low	EA
17	500-1426	Adj Collar complete - LEFT	EA	47	500-1130	Lower cross bar bracket - Low	EA
18	500-1027	Locking ring sleeve long	EA	47	500-1135	Lower cross bar bracket - Standard	EA
19	500-1041	Brake cable bracket - Complete	EA	48	500-1140	Spring for support arm - 2pc/SET	SET
20	500-1042	Brake adjustment screw - 2pc/SET	PR	49	500-1142	Cross bar bracket screw - Upper	EA
21	500-1044	Brake pad complete - LEFT	EA	50	500-1143	Cross bar bracket screw - Lower	EA
22	500-1045	Brake pad complete - RIGHT	EA	51	500-1075	Seat tube plugs - 8pc / SET	SET
23	500-1046	Brake pad spring - 2pc / SET	PR	52	500-1144	Cross bar centre screw	EA
24	500-1047	Cable Anchor - 2pc / SET	PR	--	500-2000	Wire Basket	EA
25	500-1448	Front Fork Assembly (R2)	EA	--	500-3000	Shopping bag	EA
26	500-1450	wheel 150mm - Complete	EA	--	500-4100	Cane Holder	EA
27	500-1058	Wheel bearing hardware	SET	--	500-4200	Cup Holder	EA
28	500-1059	Rear Axle Bolt	EA	--	500-4300	Tray	EA
29	500-1060	Front Axle Bolt	EA				
30	500-1062	Bottom tube plug	EA				
31	500-1070	Seat Side tube plug (set) - 4pc/SET	SET				
32	500-1201	Right Side frame - Escape Mini - Low	EA				
32	500-1200	Right Side frame - Escape Mini - Standard	EA				
33	500-1211	Left Side frame - Escape Mini - Low	EA				
33	500-1210	Left Side frame - Escape Mini - Standard	EA				

