

# PANTHERA MICRO



## THE PANTHERA MICRO

is a wheelchair for very young children. We all know that being able to get around on your own is very important, not least for children.

This is why the Panthera Micro is built to be extremely easy to drive in every possible respect.

The Panthera Micro is designed for children between one and five years of age (the youngest user we've met was 8 months old) and it's intended for indoor use.

## STANDARD FEATURES

- Chrome molybdenum chassis
- Titanium rear axle (10° camber)
- Height Adjustable Footplate
- Anti-tipper
- 16" or 18" rear wheels
- 3" Front castor wheel
- Balancing adjustable axle
- Light weight frame: starting at 5.15lbs (2.342Kg)
- No pushrims

## ACCESSORIES (Sold Separately)



Adjustable Push Handle



Side Guards



Spoke Protectors



Cushion



Brakes

## PANTHERA MICRO SPECIFICATIONS

<b>Seat width:</b>	9.4" (24 cm)	10.6" (27 cm)	9.4" (24 cm Long)	10.6" (27 cm Long)
<b>backrest height:</b>	8" (20 cm)			
<b>Seat height</b> w. 18" rear wheel (rear):	12.2" (31 cm)			
<b>Seat height</b> w. 20" rear wheel (rear):	13.4" (34 cm)			
<b>Seat depth:</b>	5.8-7.8" (15-20 cm)		7.8 - 9.8" (20-25 cm)	
<b>Seat angle:</b>	0°			
<b>Total length:</b>	24.4" (62 cm)		28.5" (72.5 cm)	
<b>Total height:</b>	20.25" (51.5c m)		21.25" (54 cm)	
<b>Total width:</b>	18" (46 cm)	19.25" (49 cm)	18" (46 cm)	19.25" (49 cm)
<b>Max User weight:</b>	66 lb (30 kg)			
<b>Total weight</b>	7.74 lb (3.52 kg)	7.84 lb (3.56 kg)	8.32 lb (3.78 kg)	8.42 lb (3.83 kg)
<b>Transfer weight</b>	5.15 lb (2.34 kg)	5.25 lb (2.39 kg)	5.44 lb (2.47 kg)	5.54 lb (2.52 kg)

**panthera**

Distributed exclusively  
in North America by



Triumph Mobility Inc.

www.triumphmobility.com | info@triumphmobility.com | 1.855.546.0711  
R\_1016\_Micro-SalesSheet\_EN | Printed in Canada

