



Why Choose Panthera Bambino 3?

Like Triumph, our partner Panthera is always thinking about ways to innovate. The Bambino 3 children's wheelchair reflects that attitude. Already an easy-to-drive, everyday chair for children with many different needs, Panthera's recent additions offer even more customization suited to children from ages 4 to 12

Here are just some of the benefits built into the Bambino 3:

Footrest flips up, moves forward and back, and adjusts vertically for shorter legs, as well as by angle.

Brake sits beside the wheel, and is easily accessed by even the most frail children from a lever inside the chair that moves out of the way for transfer

Push bar that's sturdy, quick release and easily adjustable in height and angle. The push bar is also removeable.the most frail children.

Seat width adaptors narrow the seat width by 7.5 mm on each side, growing with your child

Easily adjustable arm supports, lateral supports and side guards add customization and comfort

Head support that's cushioned and easily adjustable to the push bar and for multiple head positions for the best comfort and fit

Tray table slides into existing slots for easy use, whether for snacks or browsing on a tablet. An impressive magnet feature means you easily attach it to the wheel of the chair have it to go for the whole day!

Entire chair easily transported and fits in the trunk of a car



FEATURES

Lightweight frame: starting at 10.758 ib (4.89)

Chassis made from chrome molybdenum steel, rear axle from carbon fibre

Rear wheels from 20 to 24 inches and 5-inch front caster wheel

Stable, fixed frame for stability and adjustable balance point

OPTIONAL ACCESSORIES

Panthera Bambino 3 offers a great variety of high-quality accessories that help enhance the wheelchair comfort and performance. These can be purchased separately.



Spoke protector



Height adjustable push handle



Push bar



Side guards



Cushion

CHASSIS COLOUR



SPECIFICATIONS

Seat width	9.45" (24 cm)	10.63" (27 cm)	9.45" (24 cm)	10.63" (27 cm)	11.8" (30 cm)	13" (33 cm)
Backrest height	7.87" - 11" (20 - 28 cm) or 10.6" - 13.8" (27 - 35 cm)					
Seat height (rear)	15" (39 cm)					
Seat height (front)	16.1" (41 cm)					
Seat angle	4°					
Seat depth	7" - 9.8" (18 - 25 cm)		9.8" - 11.8" (25 - 30 cm)			
Total length	25.2" (65 cm)		25.6" - 28.1" (65 - 71.5 cm)			
Total height	22.6" - 25.6" (57.5 - 65 cm)					
Total width	21.65" (55 cm)	22.8" (58 cm)	21.65" (55 cm)		24" (61 cm)	25.2" (64 cm)
Max User weight	132 lb (60 kg)					
Total weight	21.98 lb (9.92 kg)				22.48 lb (10.22 kg)	
Transfer weight	12.26 lb (5.57 kg)		12.76 lb (5.80 kg)	13.26 lb (6.03 kg)	12.3 lb (5.6 kg)	