



TRIUMPH
escape

THE NEXT STEP IN MOBILITY

The Escape rollator offers you the independence to maintain an active lifestyle and the confidence that comes from a reliable and comfortable design.



This adjustable, foldable rollator suits almost any user, with a compact construction that makes it easy to fit into the trunk of a car. You can even fold it partially while in use to fit through narrow doorways!

Customers like the Escape because it is compact yet comfortable, with a large padded seat and ergonomic backrest. Made of lightweight aluminum, the Escape is designed for durability and its flexibility to function while walking or seated adds comfort to a user's daily

routine. Want to walk to the store and stop along the way for a chat or a coffee? No problem when you know you've always got a seat available to save your energy.

Special features include brake handles designed at a more natural angle to reduce stress on wrists, hand grips that are wider for easier grip, 3 seat heights available (19, 21, and 24 inches) and a slim profile that's easier to store or fold partially to fit through narrow passageways. A convenient, foldable shopping bag is included.

FEATURES

Ergonomic padded backrest with adjustable depth

Anatomical handles with easy-to-apply, hidden cable braking system

Large padded seat

A compact folding system activated by the lift of a release handle

Light weight, aluminum construction with strong T-frame

Large 8-inch wheels for smoother navigation



3 seat heights (19, 21, and 24 inches)

Convenient, foldable shopping bag (Shopping bag with cover available)

Slow down brake, cane holder, cup holder & tray (sold separately)

ACCESSORIES



Cane Holder



Cup Holder



Bamboo Tray



Slow down Brake

SPECIFICATIONS

Colour	Model Numbers		
Charcoal	500-10191	500-10211	500-10241
Red	500-10195	500-10215	500-10245
Blue	500-10192	500-10212	500-10242
Description	Super Low	Low	Standard
Seat height	19" (48 cm)	21" (54 cm)	24" (61 cm)
Handle height range	28 - 36" (71 - 91.5 cm)	29.5 - 37.5" (75 - 95.3 cm)	32 - 40" (81.3 - 101.6 cm)
Distance between handles	18" (45.7 cm)	18" (45.7 cm)	18" (45.7 cm)
Distance between wheels	19" (48.25 cm)	19" (48.25 cm)	19" (48.25 cm)
Seat size	18.5 x 9.5" (47 x 24 cm)	18.5 x 9.5" (47 x 24 cm)	18.5 x 9.5" (47 x 24 cm)
Overall width	25" (63.5 cm)	25" (63.5 cm)	25" (63.5 cm)
Overall length	25.6" (65 cm)	25.6" (65 cm)	25.6" (65 cm)
Depth folded	11" (28 cm)	11" (28 cm)	11" (28 cm)
Wheel diameter	8" (200 mm)	8" (200 mm)	8" (200 mm)
Weight capacity of rollator	300 lb (136 kg)	300 lb (136 kg)	300 lb (136 kg)
Weight capacity of shopping bag	11 lb (5 kg)	11 lb (5 kg)	11 lb (5 kg)
Rollator Weight	15.4 lb (7 kg)	15.6 lb (7.1 kg)	15.8 lb (7.2 kg)