



Parts List: Triumph Mini

Click [part number](#) to see part image

1. Brake Handle

	500-1510	Escape Mini Brake handle (Black) - LEFT	EA
	500-1514	Escape Mini Brake handle (Black) - RIGHT	EA
1a	500-1002	Backrest connector - LEFT	EA
1b	500-1009	Anatomical handles (1 PR)	PR
1c	500-1011	Housing / backstrap connector screw & nut	EA
1d	500-1005	Backrest connector - RIGHT	EA

2. Adj. Collar

2a	500-1025	Adj Collar complete (RED) - LEFT	EA
2a	500-1525	Adj Collar complete (GREY) - LEFT	EA
2b	500-1026	Adj Collar complete (RED) - RIGHT	EA
2b	500-1526	Adj Collar complete (GREY) - RIGHT	EA

3. Adj. Collar Sleeve *(common to both sides)*

	500-1027	Locking ring sleeve long	EA
--	--------------------------	--------------------------	----

4. T-frame

4d	500-1041	Brake cable bracket - Complete	EA
4e	500-1070	Seat Side tube plug (set) - 4pc/SET	EA
4f	500-1448	Front Fork (NEW STYLE) Assembly	EA
4g	500-1062	Bottom tube plug	EA
4h	500-1101	Basket hook - LEFT	EA
4i	500-1102	Basket hook screw	EA
4j	500-1103	Seat to frame screw	EA
4k	500-1104	Red clips for shopping bag (Left & Right)	EA
4m	500-1100	Basket hook - RIGHT	EA

5. Backstrap

5a	500-1403	Backstrap (Large Padded) - Complete	EA
5b	500-1405	Backstrap - Complete	EA
5c	500-1004	Backstrap screw & Nut - Complete - 2pc/SET	EA

6. Seat

	500-1071	Seat Mechanism - Complete (Assembled)	EA
6b	500-1085	Seat Lock Assembly	EA
6c	500-1590	Soft Seat - LEFT	EA
6d	500-1591	Soft Seat - RIGHT	EA
6e	500-1592	Soft Seat - CENTRE	EA
6f	500-1095	Seat Handle Strap	EA

7. Cross Bar

7a	500-1115	Cross bar complete - Escape Low	EA
7b	500-1130	Lower cross bar bracket - Low	EA
7c	500-1140	Spring for support arm - 2 pc/SET	EA
7d	500-1142	Cross bar bracket screw - Upper	EA
7e	500-1143	Cross bar bracket screw - Lower	EA
7f	500-1144	Cross bar screw	EA

8. Rear wheel

8a	500-1044	Brake pad complete - LEFT	EA
8b	500-1450	Wheel 150 mm - Complete	EA
8c	500-1059	Rear Axle Bolt	EA
8d	500-1045	Brake pad complete - RIGHT	EA

9. Front wheel

9a	500-1450	Wheel 150 mm - Complete	EA
9b	500-1060	Front Axle Bolt	EA

10. Shopping bag

10a	500-3000	Shopping bag - Escape	EA
	500-3100	Shopping bag with zipper - Escape	EA