

## Essentials Compact Rollator Parts List



Call out	Part #	Description	UoM
1	<a href="#">520-1000</a>	Brake handle w. cable (Right)	EA
	<a href="#">520-1001</a>	Brake handle w. cable (Left)	EA
2	<a href="#">700-3003</a>	Backrest - compact	EA
3	<a href="#">520-1005</a>	Handle lock button (Right)	EA
	<a href="#">520-1006</a>	Handle lock button (Left)	EA
4	<a href="#">700-3005</a>	Seat	EA
5	<a href="#">520-1015</a>	Bag fixation w. seat cap (Right)	EA
	<a href="#">520-1016</a>	Bag fixation w. seat cap (Left)	EA
6	<a href="#">520-1020</a>	Locked hook	EA
7	<a href="#">700-3010</a>	Brake pad (Right)	EA
	<a href="#">700-3011</a>	Brake pad (Left)	EA
8	<a href="#">700-3015</a>	Curb climber (Right)	EA
	<a href="#">700-3016</a>	Curb climber (Left)	EA
9	<a href="#">520-1040</a>	Rear wheel	EA
10	<a href="#">520-1045</a>	Front wheel and fork	EA
11	<a href="#">700-3100</a>	Cane holder (compact)	EA
12	<a href="#">520-1200</a>	Bottle holder	EA

Click [part #](#) to see parts images