

# Specifications

Colour	Carbon Fibre
Seat height	22"
Handle height range	34.5" - 39.75"
Distance between handles	16.5"
Distance between rear wheels	19.5"
Seat size	19" x 8.25"
Overall width	23.75"
Overall length	25"
Depth folded	10"
Wheel diameter (front)	8"
Wheel diameter (rear)	7"
Weight capacity of rollator	300 lbs
Weight capacity of shopping bag	11 lbs
Rollator weight	12.7 lbs
Backstrap, Cane holder, Cup holder	Included

## Accessories (Sold separately)



**Slow down brakes**  
(Set of 2)



**Phone holder**  
(Universal)




Limited lifetime warranty  
[triumphmobility.com](http://triumphmobility.com)

Authorized Dealer

## Contact us

Triumph Mobility  
151 Bentley Ave, Unit 4  
Nepean, Ontario, K2E 6T7  
Tel: +1 (855) 546 0711  
Fax: +1 (866) 225 3059  
Email: [info@triumphmobility.com](mailto:info@triumphmobility.com)  
[www.triumphmobility.com](http://www.triumphmobility.com)

      TriumphMobility



SCAN ME!

©2023 Triumph Mobility Inc. All rights reserved. Triumph™ is the trademark of Triumph Mobility.

Printed in Canada

TRIUMPH  
maverick





**The lightweight  
Maverick lets you  
go anywhere. And  
everywhere.**

With the Maverick, get the support you need.

Whether you're running errands or visiting a friend, the Maverick's lightweight design gets you moving and lets you rest when you need to. Just go!

## Features

01

### Lightweight

Created with the lightest material – carbon fibre – the Maverick is just 13.8 lbs.

02

### Compact

Designed to suit an active lifestyle, The Maverick moves along different surfaces, into the car trunk or the cargo of the plane – you decide.

03

### Safe & Comfortable

Our focus on safety includes using quality materials, and offering large 8-inch wheels and a comfortable 19 x 8.25 inch seat.

04

### Accessorized

Accessories include a shopping bag, cupholder and cane holder – all included with purchase.

05

### Stylish

Sleek black and grey colouring that blends into any environment.

06

### Adjustable

Handle heights range between 34.5 and 39.75 inches suits many heights.



## Mobility freedom for everyone

Triumph Mobility was founded on a mission: to improve the daily lives of everyone we serve, helping them to move through their world with dignity and comfort.

Customers asked us for a super light rollator that would get them moving. The carbon fibre Maverick is our answer.