

TRIUMPH essentials



Triumph Essentials Products

This collection is dedicated to enhancing the daily lives of individuals by providing reliable and practical mobility solutions.

Essentials Walker: The Essentials Walker is a dependable walking companion. With its thoughtful design, it ensures safety and confidence during daily activities.

Essentials Rollator: The Essentials Rollator is designed to provide stability and ease of movement. Whether you're navigating indoors or outdoors, this rollator offers comfort and support.

Essentials Compact Rollator: A compact version of the rollator, perfect for those who need a lightweight and portable mobility aid. It's ideal for tight spaces and travel.

Essentials Shopping Rollator: A versatile shopping rollator with padded seat and spacious shopping bags, perfect for running errands.

Triumph Mobility's commitment to quality, comfort and safety continues through the Essentials Product Line. These mobility solutions are crafted to empower individuals, allowing them to move freely and maintain their independence.



Essentials Walker



Essentials Rollator



Essentials Compact Rollator



Essentials Shopping Rollator

Essential Walker

Basic Folding Walker with seat

Portable folding travel walker with 5" wheels, provides light-weight, simplicity and functionality. Padded seat, backrest, cane holder and cup holder are included.

Features



Lightweight



Adjustable handle height and seat height



Compact folding system




5" Front wheels



Weight capacity:
300 lbs (136 kgs)



Designed in Europe



Specifications



	Essential Walker (700-1000)			Essential Rollator (700-1500)		Essential Compact Rollator (700-2000)	Essential Shopping Rollator (700-5000)
Seat height	19.5" (49.5 cm)	20.5" (52 cm)	21.5" (54.6 cm)	22" (55.88 cm)	23" (58.4 cm)	20" (51 cm)	21" (53 cm)
Handle height	33.5 - 36.5" (85-93 cm)	34.5 - 37.5" (87.5-95.25 cm)	33.5 - 38.5" (90-98 cm)	35 - 38" (89-96.5 cm)	36 - 39" (91.4-99 cm)	32 - 36.75" (81.3 - 93.4 cm)	37" - 39.75" (94 - 101 cm)
Seat size	15.5 x 13.5" (39.4 x 34.3 cm)			15.5 x 13.5" (39.4 x 34.3 cm)		17 x 8" (43.2 x 20.3 cm)	14 x 9.5" (35.6 x 24 cm)
Distance between handles	16" (40.6 cm)			16.5" (42 cm)		16" (40.6 cm)	
Distance between wheels	18.5" (47 cm)			18.25" (46.4 cm)		17" (43.2 cm)	13.75" (35 cm)
Overall width	22.5" (57 cm)			24.75" (63 cm)		24.5" (62.2 cm)	20.5" (52 cm)
Overall length	20" (50.8 cm)			22.5" (57.15 cm)		25" (63.5 cm)	28.25" (71.75 cm)
Depth folded	5" (12.7 cm)			8" (20 cm)		8" (20 cm)	12" (30.5 cm)
Wheel diameter	5" (12.7 cm)			6" (15 cm)		Front: 8" (20 cm) Rear: 7" (17.75 cm)	Front: 8" (20 cm) Rear: 7" (17.75 cm)
Shopping bag weight capacity				11 lbs (5 kgs)		11 lbs (5 kgs)	Front: 33 lbs (15 kgs) Rear: 11 lbs (5 kgs)
Weight capacity	300 lbs (136 kgs)			300 lbs (136 kgs)		300 lbs (136 kgs)	300 lbs (136 kgs)
Rollator weight	12.8 lbs (5.8 kgs)			15.25 lbs (6.9 kgs)		16.95 lbs (7.7 kgs)	18 lbs (8.2 kgs)



Contact us

Triumph Mobility

151 Bentley Ave, Unit 4, Ottawa, Ontario, K2E 6T7

Tel: +1 (855) 546 0711

Fax: +1 (866) 225 3059

Email: info@triumphmobility.com

www.triumphmobility.com

The Essentials products are on limited 2 year warranty



©2024 Triumph Mobility Inc. All rights reserved. Triumph™ is the trademark of Triumph Mobility.

Printed in Canada